

Preparing Standard Acid Solution


Hydrochloric acid and sodium hydroxide are the most common strong acids and bases used in the laboratory. Both reagents need to be standardized to learn their exact concentrations. Section 10-5 in the textbook provides background information for the procedures described below.

Reagents

Methyl Red indicator

Concentrated (37 wt%) HCl

Primary standard: sodium tetraborate $\text{Na}_2\text{B}_4\text{O}_7 \times 10 \text{ H}_2\text{O}$ (“borax”). FW 381.422; Solid colorless prismatic crystals. Specific gravity is 1.72. Borax melts when heated losing water. Non-hygroscopic. Reacts with acid:



Standardizing HCl

1. Use the table inside the cover of the textbook to calculate the volume of ~37 wt% HCl that should be added to 1 L of distilled water to produce 0.1 M HCl and prepare this solution.
2. Prepare primary standard grade sodium tetraborate.
3. Weigh four samples, each containing enough $\text{Na}_2\text{B}_4\text{O}_7 \times 10 \text{ H}_2\text{O}$ to react with ~25 mL of 0.1 M HCl (around 0.5 g) and place each in a 125-mL flask. When you are ready to titrate each one, dissolve it in ~50mL of distilled water. Add 3 drops of methyl red indicator and titrate one to an orange-red color.

Perform one blank titration of 50 mL of 0.05 M NaCl containing 3 drops of indicator. Subtract the volume of HCl needed for the blank from that required to titrate $\text{Na}_2\text{B}_4\text{O}_7 \times 10 \text{ H}_2\text{O}$.

Calculate the mean HCl molarity, standard deviation, and relative standard deviation.

Example: 27.65 mL of HCl were used to titrate 0.4916 g of Borax.

Molarity of HCl is:

$$C_{\text{HCl}} = \frac{0.4916 \times 2 \times 1000}{27.65 \times 381.42} = 0.0932_3$$

Using a pH Electrode for an Acid-Base Titration

In this experiment you will use a pH electrode to follow the course of an acid-base titration. You will observe how pH changes slowly during most of the reaction and rapidly near the equivalence point. You will compute the first and second derivatives of the titration curve to locate the end point.




Figure 1. Locating the maximum position of the first derivative of a titration curve.

Analysis of a Mixture of Carbonate and Bicarbonate

Obtain the unknown in a clean 250-mL volumetric flask. Dilute to the mark and mix well. Titrate a 25.00-mL aliquot with standard HCl solution. Use the resulting titration curve to select pH suitable for end point detection, and perform triplicate titrations with these.

Identify the solute species in the unknown, and report the mass/volume percent of each. Calculate the approximate dissociation constant that can be obtained for any carbonate-containing species from the titration data.